

SINCE 1925

KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS

1073B TYVEK® AUTOCLAVE BAG WITH STEAM INDICATORS



TYPICAL APPLICATION & CHARACTERISTICS:

Designed for Steam Sterilization of Pharmaceutical Production Equipment and Components

- 100% 1073B Dupont™ Tyvek® construction allows for optimal steam penetration during sterilization
- Excellent impact resistance
- Consistently strong factory seals
- Exceptional microbial barrier
- Resists tearing and punctures from sharp corners and heavy machine parts
- Manufactured in ISO certified cleanrooms
- Lot traceable
- Certificates of Conformance sent with all shipments
- Made in USA

MATERIAL SPECIFICATIONS:

PROPERTY	UNITS	TEST METHOD	TYPICAL VALUES 1073B TYVEK®
Seal Strength	lbs/inch/width	F-88	>4.5
Basis Weight	oz/yd ²	ASTM D3776 EN ISO 536	2.2
Total Thickness	Mils	D 1777-64	7
Microbial Barrier	Log Reduction Value (LRV)	ASTM F 1608	5.2
Bendsten Air Permeability	mL/min	ISO 5636-3	572
Moisture Vapor Transmission Rate	g/m ² /24 hr	TAPPI T523	1615
Hydrostatic Head	Inches h ₂ O	AATCC TM 127 EN20811	58
Tensile Strength MD	lb/in.	ASTM D5035 EN ISO 1924-2	44
Tensile Strength CD	lb/in.	ASTM D5035 EN ISO 1924-2	45
Elongation MD	%	ASTM D5035 EN ISO 1924-2	20
Elongation CD	%	ASTM D5035 EN ISO 1924-2	24
Elmendorf Tear, MD	lbs	ASTM 1424 EN 21974	0.7
Elmendorf Tear, CD	lbs	ASTM 1424 EN 21974	0.8
Mullen Burst	psi	ASTM D774 ISO 2758	176
Spencer Puncture	in.-lb /in 2	ASTM 3420	50
Opacity	%	TAPPI T425 ISO 2471	91

SINCE 1925

KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS

1073B TYVEK® AUTOCLAVE BAG WITH STEAM INDICATORS



KEYSTONE PART NUMBER	DIMENSIONS	KEYSTONE PART NUMBER	DIMENSIONS
COVMGTVKUB1X14	1" x 14"	COVMGTVKUB18X24	18" x 24"
COVMGTVKUB3X6	3" x 6"	COVMGTVKUB18X26	18" x 26"
COVMGTVKUB5X25	5" x 25"	COVMGTVKUB18X30	18" x 30"
COVMGTVKUB6X8	6" x 8"	COVMGTVKUB18X40	18" x 40"
COVMGTVKUB6X16	6" x 16"	COVMGTVKUB18X68	18" x 68"
COVMGTVKUB6X22	6" x 22"	COVMGTVKUB20X20	20" x 20"
COVMGTVKUB6X34	6" x 34"	COVMGTVKUB20X32	20" x 32"
COVMGTVKUB8X9	8" x 9"	COVMGTVKUB24X24	24" x 24"
COVMGTVKUB8X12	8" x 12"	COVMGTVKUB24X44	24" x 44"
COVMGTVKUB8X14	8" x 14"	COVMGTVKUB25X37	25" x 37"
COVMGTVKUB8X30	8" x 30"	COVMGTVKUB26X28	26" x 28"
COVMGTVKUB10X12	10" x 12"	COVMGTVKUB26X30	26" x 30"
COVMGTVKUB10X20	10" x 20"	COVMGTVKUB30X30	30" x 30"
COVMGTVKUB10X25	10" x 25"	COVMGTVKUB30X32	30" x 32"
COVMGTVKUB12X12	12" x 12"	COVMGTVKUB32X34	32" x 34"
COVMGTVKUB12X16	12" x 16"	COVMGTVKUB34X34	34" x 34"
COVMGTVKUB12X24	12" x 24"	COVMGTVKUB34X40	34" x 40"
COVMGTVKUB12X32	12" x 32"	COVMGTVKUB36X36	36" x 36"
COVMGTVKUB12X34	12" x 34"	COVMGTVKUB36X42	36" x 42"
COVMGTVKUB14X22	14" x 22"	COVMGTVKUB36X48	36" x 48"
COVMGTVKUB14X28	14" x 28"	COVMGTVKUB36X60	36" x 60"
COVMGTVKUB14X42	14" x 42"	COVMGTVKUB38X48	38" x 48"
COVMGTVKUB15X18	15" x 18"	COVMGTVKUB40X48	40" x 48"
COVMGTVKUB16X24	16" x 24"	COVMGTVKUB42X46	42" x 46"
COVMGTVKUB16X36	16" x 36"	COVMGTVKUB48X28	48" x 28"
COVMGTVKUB16X44	16" x 44"	COVMGTVKUB48X48	48" x 48"
COVMGTVKUB18X22	18" x 22"		