

SINCE 1925

# KEYSTONE

## SAFETY

# POLYPROPYLENE BOUFFANT CAP

### 110NWI-10

## DESCRIPTION

Polypropylene O.R. Bouffant cap is made with a light weight breathable fabric that provides economical protection for hair containment. It is the best low cost solution to protect workers against dry particles, dirt, and grime and provides excellent improved hair containment.

## CHARACTERISTICS

- Spunbond polypropylene
- 100% latex free
- Provides sanitary protection
- Cool and comfortable



**SIZES AVAILABLE:** 19", 21", 24", 26", 28", 32"

**COLORS AVAILABLE:** White, Blue, Green, Yellow, Red, Orange, Navy, Custom

**PACKAGING:** 100/bag - 10 bags/case = 1,000 pieces

**COUNTRY OF ORIGIN:** China

All Keystone Safety Product components comply with the provisions of the Federal Food, Drug and Cosmetic Act, which allow repeated use in direct contact with food. Fabric conforms to FDA code of Federal Regulations Title 21, Rule 110.10 & 177.1950

ITEM	SIZE	COUNT	WT/BOX	DIMENSION	CUBE	QUANTITY/ SKID	SKID CONFIG TIERS x ROWS
110NWI-10	19"	1,000	4.3	12x9x9	.017	112	14x8
110NWI-10	21"	1,000	5	13x11x10	.023	84	12x7
110NWI-10	24"	1,000	7	13x11x11	.025	84	12x7
110NWI-10	26"	1,000	8.8	13x13x11	.029	48	8x6
110NWI-10	28"	1,000	9	13x12x13	.031	48	8x6
110NWI-10	32"	1,000	12.5	17x15x14	.056	48	8x6

1591 HYLTON ROAD PENNSAUKEN, NJ 08110  
KEYSTONESAFETY.COM | INFO@KEYSTONESAFETY.COM

TOLL FREE (800) 663.5439  
TOLL FREE FAX (800) 663.5432  
PHONE (856) 663.5740  
FAX (856) 663.6075

SINCE 1925

# KEYSTONE

## SAFETY