

SINCE 1925

# KEYSTONE

SAFETY ⇄ CLEANROOM PRODUCTS

## Autoclave Bag with Steam Indicator, Made with 1422A Tyvek®



### Description

- Designed for steam sterilization of pharmaceutical production equipment
- Cost Effective
- Tyvek® on both sides for maximum penetration during sterilization
- Fast and consistent for equipment preparation area
- Creates and consistent for equipment preparation area
- Hydrophobic
- Resists tearing and punctures
- Superior microbial barrier
- Lint-Free
- Manufactured and packaged in ISO certified cleanroom
- Made in the USA
- Certificates of Conformance sent with shipment

Tyvek® is a trademark of affiliates of DuPont de Nemours, Inc.

SINCE 1925

# KEYSTONE

SAFETY ⇄ CLEANROOM PRODUCTS

PROPERTY	UNITS	TEST METHOD	TYPICAL VALUES (422A TYVEK®)
Basis Weight	g/m <sup>2</sup>	EN ISO 536	41
Total Thickness		EN 20534	130
Air Permeability	s	ISO 5636-5	20
Water Vapor Permeability	G/M <sup>2</sup> /24 HR	ASTM 398	1700
Hydrostatic Head	cm h <sub>2</sub> O	ISO 811	120
Trapezoidal Tear MD/CD	N	ISO 9073-4	26.1/30.6
Linting Propensity	Rating	BS 6909 ( Shirley Method 21)	Pass Rating: Excellent
Burst Resistance	kPa	ISO 2960 ( 50 cm <sup>2</sup> )	108
Puncture Resistance	N	En 863	10.8
Melting Point	C		135

Keystone Part Number	Dimensions	Keystone Part Number	Dimensions
BGT2A0205	2" x 5"	BGT2A1234	12" x 34"
BGT2A0404	4" x 4"	BGT2A1240	12" x 40"
BGT2A0409	4" x 9"	BGT2A1246	12" x 46"
BGT2A0416	4" x 16"	BGT2A1248	12" x 48"
BGT2A0820	8" x 20"	BGT2A1624	16" x 24"
BGT2A09-512	9.5" x 12"	BGT2A1830	18" x 30"
BGT2A1015	10" x 15"	BGT2A2020	20" x 20"
BGT2A1020	10" x 20"	BGT2A2024	20" x 24"
BGT2A1212	12" x 12"	BGT2A2030	20" x 30"
BGT2A1216	12" x 16"	BGT2A2040	20" x 40"
BGT2A1218	12" x 18"	BGT2A2224	22" x 24"

Keystone Part Number	Dimensions
BGT2A2424	24" x 24"
BGT2A2436	24" x 36"
BGT2A2448	24" x 48"
BGT2A2844	28" x 44"
BGT2A3036	30" x 36"
BGT2A3636	36" x 36"
BGT2A3648	36" x 48"
BGT2A3660	36" x 60"
BGT2A4848	48" x 48"