

SINCE 1925

# KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS

## BLUE BHD SELF-SEALING AUTOCLAVABLE BAG



### TYPICAL APPLICATION & CHARACTERISTICS:

- Designed for steam sterilization
- Incremental improvement compared to cutting and wrapping with Bioshield™ and autoclave tape
- Avoid sterilization re-validation
- Optional steam sterilization indicator
- Eliminates need to cut material and the use of particulate generating autoclave tape
- Fast and consistent self-seal closure
- Creates secondary level of protection during aseptic presentation and setup of equipment
- Resists tearing and punctures
- 2-ply for good microbial barrier
- Manufactured and packaged in ISO certified Cleanroom
- Certificates of Conformance sent with all shipments
- Made in USA

### MATERIAL SPECIFICATIONS:

PROPERTY	UNITS	BHD MATERIAL
Basis Weight	Gm/m2	57
Dry Tensile Strength (MD)	grams	8,100
Dry Tensile Strength (CD)	grams	7,100
Elongation Dry	%	27.3
Air Permeability	liters / 100 sq. cm	110
Suter Water Head (pressure test)	mm	400
Dry Mullen Burst	Grams / sq.cm	2,150
Flammability	seconds	8

SINCE 1925

# KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS

# BLUE BHD SELF-SEALING AUTOCLAVABLE BAG



KEYSTONE PART NUMBER	DIMENSIONS
BNW0642SSSD-2P	6" X 42"
BNW1212SSSD-2P	12" X 12"
BNW1624SSSD-2P	16" X 24"
BNW1644SSSD-2P	16" X 44"
BNW1830SSSD-2P	18" X 30"
BNW1854SSSD-2P	18" X 54"
BNW2436SSSD-2P	24" X 36"
BNW4648SSSD-2P	46" X 48"