

SINCE 1925

# KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS

## BLUE BHD STOPPER BOWL COVER



Save time by using the stopper bowl cover instead of wrapping with sheets of material and autoclave tape. When used in conjunction with a secondary autoclave bag, it creates an effective wrapping system. After the outer bag is removed, this stopper bowl cover should remain in place, providing protection from microbial and particulate contamination while the remaining set-up activities are performed.

### TYPICAL APPLICATION & CHARACTERISTICS:

- Steam Sterilization Indicator
- Autoclavable up to 135°C
- Incremental improvement compared to cutting and wrapping with Bioshield™ and autoclave tape
- Eliminates need to cut material and the use of particulate generating autoclave tape
- Fast and consistent
- Creates secondary level of protection during aseptic presentation and setup of equipment
- Avoid sterilization Re-validation
- Manufactured and packaged in ISO certified Cleanroom
- Conveniently packaged: 1 per heat-sealed labeled poly bag
- Certificates of Conformance sent with all shipments
- Made in USA

SINCE 1925

# KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS

# BLUE BHD STOPPER BOWL COVER



## MATERIAL SPECIFICATIONS:

PROPERTY	UNITS	BHD MATERIAL
Basis Weight	Gm/m2	57
Dry Tensile Strength (MD)	grams	8,100
Dry Tensile Strength (CD)	grams	7,100
Elongation Dry	%	27.3
Air Permeability	liters / 100 sq. cm	110
Suter Water Head (pressure test)	mm	400
Dry Mullen Burst	Grams / sq.cm	2,150
Flammability	seconds	8

KEYSTONE PART NUMBER	DIMENSIONS
BWLNBW20	Fits 11" to 16" Diameter
BWLNBW24	Fits 12" to 20" Diameter
BWLNBW26	Fits 14" to 22" Diameter
BWLNBW30	Fits 20" to 26" Diameter
BWLNBW36	Fits 26" to 30" Diameter
BWLNBW42	Fits 28" to 35" Diameter
BWLNBW46	Fits 30" to 38" Diameter
BWLNBW48	Fits 30" to 40" Diameter
BWLNBW54	Fits 35" to 45" Diameter
BWLNBW60	Fits 40" to 50" Diameter