

SINCE 1925

KEYSTONE

SAFETY ⇄ CLEANROOM PRODUCTS

Bowl Cover, Made with Gray Tyvek®

UV Light Blocking Tyvek® FC Bowl Covers are designed to keep cleaned pharmaceutical equipment and supplies clean and protected from external sources of contamination while also providing in-process protection to light sensitive biotech products. Tyvek® is a trademark of affiliates of DuPont de Nemours, Inc.



Typical Application & Characteristics:

- Each cover is individually poly bagged and labeled for easy identification
- UV Light Blocking Tyvek® FC fabric protects light sensitive Biotech products
- Hydrophobic
- Resists tearing and punctures
- Anti-static
- Manufactured and packaged in ISO certified cleanrooms
- Certificates of Conformance sent with all shipments
- Made in USA

PROPERTY	UNITS	TEST METHOD	TYPICAL VALUES TYVEK® 400 FC
Basis Weight	oz/yd ²	ASTM D3776	1.6 oz/yd ²
Thickness	mils	ASTM D1117	7.4 mils
Breaking Strength-Grab (MD)	f	ASTM D5034	20 lbf
Breaking Strength-Grab (CD)	f	ASTM D5034	26 lbf
Hydrostatic Head	in H ₂ O	AATCC 127	4.2 in H ₂ O
Tear Resistance MD	f	ASTM D5733	5lbf
Tear Resistance CD	f	ASTM D5733	7lbf
Flammability		16 CFR 1610	Class 1

SINCE 1925

KEYSTONE

SAFETY ⇄ CLEANROOM PRODUCTS



DIMENSIONS	KEYSTONE PART NO.
Fits 11" to 16" Diameter	BWLTG2A20
Fits 12" to 20" Diameter	BWLTG2A24
Fits 14" to 22" Diameter	BWLTG2A26
Fits 20" to 26" Diameter	BWLTG2A30
Fits 26" to 30" Diameter	BWLTG2A36
Fits 28" to 35" Diameter	BWLTG2A42
Fits 30" to 38" Diameter	BWLTG2A46
Fits 30" to 40" Diameter	BWLTG2A48
Fits 35" to 45" Diameter	BWLTG2A54
Fits 40" to 50" Diameter	BWLTG2A60
Fits 50" to 60" Diameter	BWLTG2A72