

SINCE 1925

KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS

Drawstring Closure Autoclave Bag with Steam Indicator, Made with 1422A Tyvek®

Description

- Designed for steam sterilization
- Tyvek® on both sides for maximum breathability during sterilization
- Eliminates need to cut material and the use of autoclave tape
- Fast and consistent for equipment preparation area
- Creates secondary level of protection during aseptic presentation of equipment
- Sealed seams eliminate sewing threads and greatly reduce particulates
- Resists tearing and punctures
- Superior microbial barrier
- Manufactured and packaged in ISO certified cleanrooms
- Certificates of Conformance sent with all shipments
- Made in USA



SINCE 1925

KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS

SINCE 1925



SAFETY ↔ CLEANROOM PRODUCTS

PROPERTY	UNITS	TEST METHOD	TYPICAL VALUES (M22A TYVEK®)
Basis Weight	g/m ²	EN ISO 536	41
Total Thickness		EN 20534	130
Air Permeability	s	ISO 5636-5	20
Water Vapor Permeability	G/M ² /24 HR	ASTM 398	1700
Hydrostatic Head	cm h ₂ O	ISO 811	120
Trapezoidal Tear MD/CD	N	ISO 9073-4	26.1/30.6
Linting Propensity	Rating	BS 6909 (Shirley Method 21)	Pass Rating: Excellent
Burst Resistance	kPa	ISO 2960 (50 cm ²)	108
Puncture Resistance	N	En 863	10.8
Melting Point	C		135

Keystone Part Number	Size W x L	Keystone Part Number	Size W x L
BGT2A0101D	1" X 1"	BGT2A0203D	2" X 3"
BGT2A0105D	1" X 5"	BGT2A0205D	2" X 5"
BGT2A0106D	1" X 6"	BGT2A0208D	2" X 8"
BGT2A0109D	1" X 9"	BGT2A0211D	2" X 11"
BGT2A0111D	1" X 11"	BGT2A0304D	3" X 4"
BGT2A0117D	1" X 17"	BGT2A0305D	3" X 5"
BGT2A0120	1" X 20"	BGT2A0306D	3" X 6"
BGT2A0124D	1" X 24"	BGT2A0308D	3" X 8"
BGT2A0202D	2" X 2"	BGT2A0309D	3" X 9"

SINCE 1925



SAFETY ↔ CLEANROOM PRODUCTS

Keystone Part Number	Size W x L	Keystone Part Number	Size W x L
BGT2A0310D	3" X 10"	BGT2A0643D	6" X 43"
BGT2A0404D	4" X 4"	BGT2A0662D	6" X 62"
BGT2A0407D	4" X 7"	BGT2A0710D	7" X 10"
BGT2A0409D	4" X 9"	BGT2A0717D	7" X 17"
BGT2A0412D	4" X 12"	BGT2A0727D	7" X 27"
BGT2A0425D	4" X 25"	BGT2A0808D	8" X 8"
BGT2A0430D	4" X 30"	BGT2A0809D	8" X 9"
BGTA2A0443D	4" X 43"	BGT2A0812D	8" X 12"
BGT2A0448D	4" X 48"	BGT2A0817D	8" X 17"
BGT2A0606D	6" X 6"	BGT2A0824D	8" X 24"
BGT2A0608	6" X 8"	BGT2A0831D	8" X 31"
BGT2A0612D	6" X 12"	BGT2A0843D	8" X 43"
BGTA2A0617D	6" X 17"	BGT2A0856D	8" X 56"
BGT2A0623D	6" X 23"	BGT2A0865D	8" X 65"
BGT2A0633	6" X 33"	BGT2A0872D	8" X 72"
BGT2A0638D	6" X 38"	BGT2A09512D	9.5" X 12"

Keystone Part Number	Size W x L	Keystone Part Number	Size W x L
BGT2A1015D	10" X 15"	BGT2A1248D	12" X 48"
BGT2A1020D	10" X 20"	BGT2A1416D	14" X 16"
BGT2A1024D	10" X 24"	BGT2A1418D	14" X 18"
BGT2A1027D	10" X 27"	BGT2A1420D	14" X 20"
BGT2A1030D	10" X 30"	BGT2A1422D	14" X 22"
BGT2A1037D	10" X 37"	BGT2A1431D	14" X 31"
BGT2A1043D	10" X 43"	BGT2A1519D	15" X 19"
BGT2A1057D	10" X 57"	BGT2A155195D	15.5" X 19.5"
BGT2A1062D	10" X 62"	BGT2A1526D	15" X 26"
BGT2A1212D	12" X 12"	BGT2A1530D	15" X 30"
BGT2A1216D	12" X 16"	BGT2A1544D	15" X 44"
BGT2A1218D	12" X 18"	BGT2A1616D	16" X 16"
BGT2A1220D	12" X 20"	BGT2A1618D	16" X 18"
BGT2A1224D	12" X 24"	BGT2A1624D	16" X 24"
BGT2A1227D	12" X 27"	BGT2A1627D	16" X 27"
BGT2A1236D	12" X 36"	BGT2A1630D	16" X 30"

Keystone Part Number	Size W x L	Keystone Part Number	Size W x L
BGT2A1632D	16" X 32"	BGT2A2033D	20" X 33"
BGT2A1638D	16" X 38"	BGT2A2037D	20" X 37"
BGT2A1642D	16" X 42"	BGT2A2040D	20" X 40"
BGT2A1819D	18" X 19"	BGT2A2218D	22" X 18"
BGT2A1824D	18" X 24"	BGT2A2222D	22" X 22"
BGT2A1830D	18" X 30"	BGT2A2225D	22" X 25"
BGT2A1836D	18" X 36"	BGT2A2230D	22" X 30"
BGT2A1848D	18" X 48"	BGT2A2325D	23" X 25"
BGT2A1854D	18" X 54"	BGT2A2424D	24" X 24"
BGT2A1855D	18" X 55"	BGT2A2430D	24" X 30"
BGT2A1870D	18" X 70"	BGT2A2436D	24" X 36"
BGT2A2020D	20" X 20"	BGT2A2440D	24" X 40"
BGT2A2022D	20" X 22"	BGT2A2442D	24" X 42"
BGT2A2024D	20" X 24"	BGT2A2448D	24" X 48"
BGT2A2026D	20" X 26"	BGT2A2455D	24" X 55"
BGT2A2030D	20" X 30"	BGT2A2462D	24" X 62"

Keystone Part Number	Size W x L	Keystone Part Number	Size W x L
BGT2A2632D	26" X 32"	BGT2A3258D	32" X 58"
BGT2A2634D	26" X 34"	BGT2A3263D	32" X 63"
BGT2A2642D	26" X 42"	BGT2A3434D	34" X 34"
BGT2A2818D	28" X 18"	BGT2A3437D	34" X 37"
BGT2A2828D	28" X 28"	BGT2A3441D	34" X 41"
BGT2A2834D	28" X 34"	BGT2A3448D	34" X 48"
BGT2A2838D	28" X 38"	BGT2A3636D	36" X 36"
BGT2A2844D	28" X 44"	BGT2A3638D	36" X 38"
BGT2A3030D	30" X 30"	BGT2A3643D	36" X 43"
BGT2A3036D	30" X 36"	BGT2A36456D	36" X 46"
BGT2A3040D	30" X 40"	BGT2A3648D	36" X 48"
BGT2A3045D	30" X 45"	BGT2A3664D	36" X 64"
BGTA2A3053D	30" X 53"	BGT2A4040D	40" X 40"
BGT2A3060D	30" X 60"	BGT2A4047D	40" X 47"
BGT2A3246D	32" X 46"	BGT2A4049D	40" X 49"
BGT2A3253D	32" X 53"	BGT2A4052D	40" X 52"

SINCE 1925



SAFETY ⇄ CLEANROOM PRODUCTS

Keystone Part Number	Size W x L	Keystone Part Number	Size W x L
BGT2A4545D	45" X 45"	BGT2A4860D	48" X 60"
BGT2A4848D	48" X 48"	BGT2A4864D	48" X 64"
BGT2A4850D	48" X 50"		

SINCE 1925



SAFETY ⇄ CLEANROOM PRODUCTS