

SINCE 1925

KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS

GREEN GHD MUFFS



Used as an autoclavable cover these Green Heavy Duty (GHD) Muffs eliminate the need for time consuming wrapping with sheets of material and autoclave tape.

TYPICAL APPLICATION & CHARACTERISTICS:

- Autoclavable up to 135°C
- Certificates of Conformance sent with shipment
- Conveniently packaged: 100 per heat-sealed, labeled poly bag
- Manufactured and packaged in ISO certified Cleanroom
- Color: Green
- Made in USA

MATERIAL SPECIFICATIONS:

PROPERTY	UNITS	GHD MATERIAL
Basis Weight	Gm/m ²	57
Dry Tensile Strength (MD)	grams	8,100
Dry Tensile Strength (CD)	grams	7,100
Elongation Dry	%	27.3
Air Permeability	liters / 100 sq. cm	110
Suter Water Head (pressure test)	mm	400
Dry Mullen Burst	Grams / sq.cm	2,150
Flammability	seconds	8

SINCE 1925

KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS

GREEN GHD MUFFS



KEYSTONE PART NUMBER	DIMENSIONS
MFGHD04	4" used for openings 1" to 1 3/4"
MFGHD05	5" used for openings 1 3/4" to 2 1/2"
MFGHD07	7" used for openings 2 1/2" to 4"
MFGHD08	8" used for openings 3 1/2" to 4 3/4"
MFGHD09	9" used for openings 4" to 5 1/2"
MFGHD12	12" used for openings 5 1/2" to 8"
MFGHD16	16" used for openings 8" to 12"
MFGHD18	18" used for openings 10" to 14"