

SINCE 1925

KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS

Heat- Sealed Steam Sterilization Pouch with Chevron Peel, Made with Tyvek®

Description



- Tyvek®/PE Film construction for excellent strength and exceptional microbial barrier
- Designed for steam sterilization
- Fiber free peelable seals
- Excellent clarity for product visibility
- Recommended for sterilization of pharmaceutical equipment
- Manufactured and packaged in ISO certified cleanroom
- Works with both impulse and continuous band sealers
- Improves the efficiency of the parts preparation department
- Certificates of Conformance sent with each shipment
- Made in USA

Tyvek® is a trademark of affiliates of DuPont de Nemours, Inc.

PROPERTY	UNITS	TEST METHOD	TYPICAL VALUES 1073B TYVEK®	TYPICAL VALUES PET/POLYOLEFIN FILM
Seal Strength	Lbs/inchwidth	F-88	>1.5	>1.5
Basis Weight	oz/yd2	ASTM D3776 EN ISO 536	2.2	N/A
Total Thickness	Ms	D 1777-64	7	2.6
Microbial Barrier	Log Reduction Value (LRV)	ASTM F 1608	5.2	N/A
Bendsten Air Permeability	mL/min	ISO 5636-3	572	N/A
Moisture Vapor Transmission Rate	g/m2/24 hr	TAPPI T523	1615	N/A
Hydrostatic Head	Inches H2O	AATCC TM 127 EN20811	58	N/A
Tensile Strength MD	lb/in.	ASTM D5035 EN ISO 1924-2	44	N/A
Tensile Strength CD	lb/in.	ASTM D5035 EN ISO 1924-2	45	N/A
Elongation MD	%	ASTM D5035 EN ISO 1924-2	20	N/A
Elongation CD	%	ASTM D5035 EN ISO 1924-2	24	N/A
Elmendorf Tear, MD	lbs	ASTM 1424 EN 21974	0.7	N/A
Elmendorf Tear, CD	lbs	ASTM 1424 EN 21974	0.8	N/A
Mullen Burst	psi	ASTM D774 ISO 2758	176	N/A
Spencer Puncture	in.-lb//in.2	ASTM 3420	50	N/A
Opacity	%	TAPPI T425 ISO 2471	91	6

Size W x L	Keystone Part No.
3" x 6"	BBMGTVKHS0306
6" x 8"	BBMGTVKHS0608
6" x 10"	BBMGTVKHS0610
6" x 12"	BBMGTVKHS0612
8" x 10"	BBMGTVKHS0810
8" x 15"	BBMGTVKHS0815
10" x 20"	BBMGTVKHS1020
10" x 30"	BBMGTVKHS1030
12" x 12"	BBMGTVKHS1212
12" x 17"	BBMGTVKHS1217
12" x 18"	BBMGTVKHS1218
12" x 24"	BBMGTVKHS1224
16" x 17"	BBMGTVKHS1617
18" x 22"	BBMGTVKHS1822
24" x 18"	BBMGTVKHS2418
24" x 24"	BBMGTVKHS2424
24" x 30"	BBMGTVKHS2430
24" x 40"	BBMGTVKHS2440
24" x 48"	BBMGTVKHS2448

