

SINCE 1925



POLYPROPYLENE ISOLATION GOWN

DESCRIPTION

Polypropylene isolation gown that is a lightweight breathable comfortable fabric that provides economical protection. It is the best low cost solution to protect workers against dry particulates, dirt and grime. Commonly used in hospital, healthcare, and medical facilities.

CHARACTERISTICS

- Non woven 100% spunbond polypropylene
- 100% latex free

SIZES AVAILABLE: One size fits all

COLORS AVAILABLE: Blue and Yellow

PACKAGING: 50/case

COUNTRY OF ORIGIN: China

ISO-NW



All Keystone Safety Product components comply with the provisions of the Federal Food, Drug and Cosmetic Act, which allow repeated use in direct contact with food. Fabric conforms to FDA code of Federal Regulations Title 21, Rule 110.10 & 177.1950

ITEM	SIZE	COUNT	WT/BOX	DIMENSION	CUBE	QUANTITY/ SKID	SKID CONFIG TIERS x ROWS
ISO-NW	OSFA	50	7.2	11x10x16	.028	60	12x5

1591 HYLTON ROAD PENNSAUKEN, NJ 08110
KEYSTONESAFETY.COM | INFO@KEYSTONESAFETY.COM

TOLL FREE (800) 663.5439
TOLL FREE FAX (800) 663.5432
PHONE (856) 663.5740
FAX (856) 663.6075

SINCE 1925

