

SINCE 1925



LAMINATED ISOLATION GOWN

DESCRIPTION

Polyethylene laminated polypropylene material that has the added benefit of a solid sheet of polyethylene bonded to the surface providing protection where liquid splash is a concern.

CHARACTERISTICS

- Excellent protection against fluids and splashes
- 100% latex free
- Extremely durable
- The polyethylene film is impervious and non-linting with great tensile strength

ISO-NWP



SIZES AVAILABLE: One size fits all

COLORS AVAILABLE: White

PACKAGING: 50/case

COUNTRY OF ORIGIN: China

All Keystone Safety Product components comply with the provisions of the Federal Food, Drug and Cosmetic Act, which allow repeated use in direct contact with food. Fabric conforms to FDA code of Federal Regulations Title 21, Rule 110.10 & 177.1950

ITEM	SIZE	COUNT	WT/BOX	DIMENSION	CUBE	QUANTITY/ SKID	SKID CONFIG TIERS x ROWS
ISO-NWP	OSFA	50	10.2	13x11x18	.041	60	12x5

1591 HYLTON ROAD PENNSAUKEN, NJ 08110
KEYSTONESAFETY.COM | INFO@KEYSTONESAFETY.COM

TOLL FREE (800) 663.5439
TOLL FREE FAX (800) 663.5432
PHONE (856) 663.5740
FAX (856) 663.6075

SINCE 1925

