


SINCE 1925

KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS

WHITE POLYPROPYLENE BOWL COVERS



The use of Keystone Cleanroom Product's Polypropylene Bowl Covers is an efficient and consistent means of protecting critical surfaces from external sources of particulate contamination after cleaning. The breathable material allows residual moisture from the cleaning process to evaporate.

- Conveniently covers cleaned pharmaceutical production equipment
- Breathable
- Manufactured and packaged in a Certified ISO8 Certified
- Latex-and talc free rubber elastic
- Low Lint
- Packaged 1 per labeled poly bag with Date of Manufacture and Lot Number for Complete Traceability

PHYSICAL PROPERTIES

Material –100% Non-woven polypropylene Fabric,
CAS# 9003-07-4

Color White
Basis Weight 1.25 oz. sq. yd

STRENGTH PROPERTIES – (MD machine direction; CD cross direction.)

Elongation (MD) 69.2%
Elongation (CD) 68.4%
Tensile (MD) 99.0 N / 5cm
Tensile (CD) 60.7 N / 5cm

THERMAL PROPERTIES

Melting point is > 225°F.

BARRIER PROPERTIES


Vapor Pressure negligible mm HG/70F

SINCE 1925

KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS

WHITE POLYPROPYLENE BOWL COVERS



KEYSTONE PART NUMBER	DIMENSIONS
BWLPPW20E	Fits 11" to 16" Diameter
BWLPPW24E	Fits 12" to 20" Diameter
BWLPPW26E	Fits 14" to 22" Diameter
BWLPPW30E	Fits 20" to 26" Diameter
BWLPPW36E	Fits 26" to 30" Diameter
BWLPPW42E	Fits 28" to 35" Diameter
BWLPPW46E	Fits 30" to 38" Diameter
BWLPPW48E	Fits 30" to 40" Diameter
BWLPPW54E	Fits 35" to 45" Diameter
BWLPPW60E	Fits 40" to 50" Diameter
BWLPPW72E	Fits 50" to 60" Diameter
BWLPPW90E	Fits 60" to 76" Diameter