

SINCE 1925

KEYSTONE

SAFETY ⇄ CLEANROOM PRODUCTS

PRE-STERILIZED CLEANROOM TRASH BAGS



Simplify the removal of waste from your aseptic cleanroom without introducing external sources of microbial contamination, by using pre-sterilized trash liners.

TYPICAL APPLICATION & CHARACTERISTICS:

- Packaged in an ISO6 certified Cleanroom
- Strong, 2 mil clear polyethylene
- Standard sizes are 24" x 28" and 36" x 55"
- Validated Sterile to a 10^{-6} SAL
- Certificates of Conformance and Sterility sent with each shipment
- Conveniently packaged in Tyvek®/Mylar chevron peel pouches: 1, 2, 3, 4, 5 or 10 per pack configurations
- Made in USA

MATERIAL SPECIFICATIONS:

PHYSICAL PROPERTIES	UNITS
Material--Low Density Polyethylene	2 mil
Dimensions	(+/- 1")
S-2428CP	24" X 28"
S-3655CP	36" X 55"
Coefficient of Friction	0.15
STRENGTH PROPERTIES	UNITS
Tensile Break	3,900 lbs /sq Inch (MD); 2,500 lbs /sq inch (CD)
1% Seacant Modulus	22,000 lbs / sq inch (MD); 28,000 lbs / sq inch (CD)
Elongation	300% (MD); 700% (CD)
Dart Impact	100 grams/ mil F50

SINCE 1925

KEYSTONE

SAFETY ⇄ CLEANROOM PRODUCTS

PRE-STERILIZED CLEANROOM TRASH BAGS



KEYSTONE PART NUMBER	DIMENSIONS
S-2428CP-1	24" X 28", 1 per pack
S-2428CP-2	24" X 28", 2 per pack
S-2428CP-5	24" X 28", 5 per pack
S-2428CP-10	24" X 28", 10 per pack
S-3655CP-1	36" X 55", 1 per pack
S-3655CP-2	36" X 55", 2 per pack
S-3655CP-3	36" X 55", 3 per pack
S-3655CP-4	36" X 55", 4 per pack
S-3655CP-5	36" X 55", 5 per pack
S-3655CP-10	36" X 55", 10 per pack