

SINCE 1925

KEYSTONE

SAFETY ⇄ CLEANROOM PRODUCTS

PRE-STERILIZED ZIP LOCK BAGS



Keystone Cleanroom Products pre-sterilized zip lock bags are designed to reduce the risks of producing false positive media plates. This system provides complete protection from external sources of microbial contamination from the time the plates are removed from the isolator all the way through incubation.

The clear film material offers excellent clarity combined with strength and durability. These bags feature an integrated slide lock zipper seal for consistent and efficient sealing using the isolator gloves.

TYPICAL APPLICATION & CHARACTERISTICS:


- Packaged in an ISO6 certified Cleanroom
- Strong, 2.0, 2.7 or 4.0 mil clear Low density polyethylene
- Validated Sterile to a 10^{-6} SAL
- Certificates of Conformance and Sterility sent with each shipment
- Single or double pouched in Tyvek®/Mylar chevron peel pouches with an additional sealed poly bag layer
- Tyvek® peel pouch packages have optional integrated hang holes for suspension of the sterile zip lock packs in the isolator during VHP decontamination
- 5, 10 or 20 per sterile pack configurations
- Made in USA

SINCE 1925

KEYSTONE

SAFETY ⇄ CLEANROOM PRODUCTS

PRE-STERILIZED ZIP LOCK BAGS



KEYSTONE PART NUMBER	DESCRIPTION
0606ZLCRGST20	6" x 6", 20 per pack, 4 mil
0808ZLCRGST20	8" x 8", 20 per pack, 4 mil
1212ZLCRGST20	12" x 12", 20 per pack, 4 mil
1215SZLDPST10	312" x 15", 10 per pack, 2.7 mil
1215SZLDPST20	12" x 15", 20 per pack, 2.7 mil
1616ZLCRGST10	16" x 16", 10 per pack, 4 mil