

SINCE 1925

KEYSTONE

SAFETY ⇄ CLEANROOM PRODUCTS

PURECLEAN™

SELF-SEALING 1421B TYVEK®

AUTOClave BAG WITH CLEAR FILM SIDE

AND STEAM INDICATOR

PureClean™ Self-Sealing Autoclave Bags are made from a transparent film and 1421B DuPont™ Tyvek®, which is a new grade of uncoated, untreated, flexible Tyvek® material specifically designed for pharmaceutical use in steam sterilization applications.



- Uncoated Tyvek® eliminates potential of residue deposits on equipment, post sterilization
- Clear Side enables visual confirmation of dryness for sterilized parts and equipment, which is required by EU GMP Annex 1
- Convenient Self-Seal Closure eliminates the need for autoclave tape or the validation and use of a heat sealer
- Fast and consistent for equipment preparation area
- Resists tearing and punctures
- Superior microbial barrier
- Manufactured and packaged in a certified ISO 6 cleanroom
- Made in the USA
- Certificates of Conformance sent with shipment

MATERIAL SPECIFICATIONS:

Property	Comparable standard test method	Units	Nominal value (range) Tyvek® 1421B
Preliminary product specifications (English units)			
Basis weight	ASTM D3776 ¹	oz/yd ²	1.22 [1.16–1.26]
Delamination	ASTMD2724 ²	lb _f	0.08 [>0.03]
Gurley Hill porosity	TAPPI T460 ³	sec/100 cc	20 [8–36]
Preliminary product specifications (metric units)			
Basis weight	ASTM D3776 ¹	g/m ²	41 [39–43]
Delamination	ASTMD2724 ²	N	0.36 [>0.13]
Gurley Hill porosity	TAPPI T460 ³	sec/100 cc	20 [8–36]

SINCE 1925

KEYSTONE

SAFETY ⇄ CLEANROOM PRODUCTS

SELF-SEALING 1421B TYVEK® AUTOCLAVE BAG WITH CLEAR FILM SIDE AND STEAM INDICATOR

Property (Continued)	Comparable standard test method	Units	Typical value Tyvek® 1421B
----------------------	---------------------------------	-------	----------------------------

Preliminary miscellaneous properties (English units)

Thickness	EN ISO 534 ¹	mils	6 (0.15mm)
Tensile strength, MD	EN ISO 1924-2 ²	lb /in	7.6
Tensile strength, CD	EN ISO 1924-2 ²	lb /in	6.4
Elmendorf tear, MD	ASTM D1424	lb _f	2.4
Elmendorf tear, CD	ASTM D1424	lb _f	1.9
Hydrostatic head	AATCC TM 127 ³	in H ₂ O	44
Elongation, MD	EN ISO 1924-2 ²	%	10
Elongation, CD	EN ISO 1924-2 ²	%	15
Puncture resistance	EN 863 ⁴	lb _f	2.7
Mullen burst	ISO 2758	psi	53
Bendsten air permeability	ISO 5636-3 ⁵	mL/min	741
Microbial barrier	ASTM F1608	LRV	~4

MD = machine direction; CD = cross direction

1. Surface 2 cm², pressure 50 kPa.

2. Modified for speed, sample width 2.54 cm (1 inch).

3. Rate of use: 60 cm H₂O/min.

4. Modified for conditioning.

5. Electronic device.

Property	Comparable standard test method	Units	Typical value Tyvek® 1421B
----------	---------------------------------	-------	----------------------------

Preliminary miscellaneous properties (metric units)

Thickness	EN ISO 534 ¹	µm	152
Tensile strength, MD	EN ISO 1924-2 ²	N	34
Tensile strength, CD	EN ISO 1924-2 ²	N	28
Elmendorf tear, MD	ASTMD1424	N	10.6
Elmendorf tear, CD	ASTMD1424	N	8.5
Hydrostatic head	AATCC TM 127 ³	cm H ₂ O	112
Elongation, MD	EN ISO 1924-2 ²	%	10
Elongation, CD	EN ISO 1924-2 ²	%	15
Puncture resistance	EN 863 ⁴	N	12
Mullen burst	ISO 2758	kPa	365
Bendsten air permeability	ISO 5636-3 ⁵	mL/min	741
Microbial barrier	ASTM F1608	LRV	~4

MD = machine direction; CD = cross direction

1. Surface 2 cm², pressure 50 kPa.

2. Modified for speed, sample width 2.54 cm (1 inch).

3. Rate of use: 60 cm H₂O/min.

4. Modified for conditioning.

5. Electronic device.

SELF-SEALING 1421B TYVEK® AUTOCLAVE BAG WITH CLEAR FILM SIDE AND STEAM INDICATOR



KEYSTONE PART NUMBER	DIMENSIONS	KEYSTONE PART NUMBER	DIMENSIONS
6STKRP0205SS	2" x 5"	6STKRP1830SS	18" X 30"
6STKRP0404SS	4" X 4"	6STKRP2020SS	20" X 20"
6STKRP0409SS	4" x 9"	6STKRP2024SS	20" x 24"
6STKRP0416SS	4" x 16"	6STKRP2030SS	20" x 30"
6STKRP0820SS	8" x 20"	6STKRP2040SS	20" x 40"
6STKRP09-512SS	9.5" x 12"	6STKRP2224SS	22" x 24"
6STKRP1015SS	10" x 15"	6STKRP2424SS	24" x 24"
6STKRP1020SS	10" X 20"	6STKRP2436SS	24" x 36"
6STKRP1212SS	12" x 12"	6STKRP2448SS	24" x 48"
6STKRP1216SS	12" x 16"	6STKRP2844SS	28" x 44"
6STKRP1218SS	12" x 18"	6STKRP3036SS	30" X 36"
6STKRP1234SS	12" x 34"	6STKRP3636SS	36" x 36"
6STKRP1240SS	12" x 40"	6STKRP3648SS	36" X 48"
6STKRP1246SS	12" x 46"	6STKRP3660SS	36" x 60"
6STKRP1248SS	12" x 48"	6STKRP4848SS	48" x 48"
6STKRP1624SS	16" x 24"	6STKRP4860SS	48" x 60"