

SINCE 1925

KEYSTONE


SAFETY ⇄ CLEANROOM PRODUCTS

PureClean™ 1421B Tyvek® Autoclavable Stopper Bowl Covers With Tunnelized Elastic

The use of Keystone Cleanroom Products PureClean™ line of Tyvek® Autoclavable Stopper Bowl Covers is an efficient and consistent means of preparing pharmaceutical production equipment for steam sterilization. PureClean™ autoclavable bowl covers do not have a water soluble antistat agent applied during its manufacturing process like that of Tyvek 1422A. The covers also are manufactured using a tunnelized sewing technique to encapsulate the latex-free, talc-free elastic. After sterilization, the covered critical surfaces are protected from particulate and microbial contamination during the aseptic installation and setup activities.

Description

- Uncoated Tyvek® for maximum breathability during sterilization without the chance of residue deposits on equipment, post sterilization
- Individually packaged in sealed, labeled poly bags for traceability and protection from contamination until the time of use
- The tunnelized elastic reduces the risk of fiber shedding
- Breathable, while providing a superior barrier post sterilization
- Tyvek® has the highest microbial barrier properties among all steam sterilization wrapping materials
- Manufactured and packaged in a ISO Certified Cleanroom
- Eliminates Autoclave Tape Residue
- Latex and Talc free rubber elastic
- Lint-Free



SINCE 1925

KEYSTONE

SAFETY ✧ CLEANROOM PRODUCTS

Fits Openings with Diameters	Keystone Part No.
12" to 20"	BCTK24T
14" to 22"	BCTK26T
20" to 26"	BCTK30T
26" to 30"	BCTK36T
28" to 23""	BCTK38T
28" to 35"	BCTK42T
30" to 40"	BCTK48T
35" to 45"	BCTK54T
40" to 50"	BCTK60T

SINCE 1925

KEYSTONE

SAFETY ✧ CLEANROOM PRODUCTS

Contact Us:

+1 (856) 663-5740

info@keystonecrp.com