

SINCE 1925

KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS

BLUE SUPER DUTY SMS STOPPER BOWL COVER



Save time by using the stopper bowl cover instead of wrapping with sheets of SMS material and autoclave tape. When used in conjunction with a secondary autoclave bag, it creates an effective wrapping system. After the outer bag is removed, this stopper bowl cover should remain in place, providing protection from microbial and particulate contamination while the remaining set-up activities are performed.

TYPICAL APPLICATION & CHARACTERISTICS:

- Eliminates autoclave tape residue
- SMS has excellent microbial barrier properties
- Can be sterilized in **135°** autoclave cycles
- Creates secondary level of protection during aseptic presentation and setup of equipment
- Hydrophobic
- Resists tearing and punctures
- Lint-free
- Manufactured and packaged in ISO certified cleanrooms
- Conveniently packaged: 1 per heat-sealed and labeled poly bag
- Certificates of Conformance sent with all shipments
- Made in USA

MATERIAL SPECIFICATIONS:


| PROPERTY | SPECIFICATION | UNIT |
|-------------------------------|---------------|-------|
| Basis Weight | 88 | Gm/m2 |
| Grab Strength (CD) | 40 | lbs |
| Tear Strength (MD) | 11 | lbs |
| Mullen Burst Strength (Gross) | 60 | psi |
| Air Permeability | 25 | cfm |
| Hydrohead | 70 | cm |

SINCE 1925

KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS

BLUE SUPER DUTY SMS STOPPER BOWL COVER



| KEYSTONE PART NUMBER | DIMENSIONS |
|----------------------|--------------------------|
| BWL-KCRP-600-20 | Fits 11" to 16" Diameter |
| BWL-KCRP-600-24 | Fits 12" to 20" Diameter |
| BWL-KCRP-600-26 | Fits 14" to 22" Diameter |
| BWL-KCRP-600-30 | Fits 20" to 26" Diameter |
| BWL-KCRP-600-36 | Fits 26" to 30" Diameter |
| BWL-KCRP-600-42 | Fits 28" to 35" Diameter |
| BWL-KCRP-600-46 | Fits 30" to 38" Diameter |
| BWL-KCRP-600-48 | Fits 30" to 40" Diameter |
| BWL-KCRP-600-54 | Fits 35" to 45" Diameter |
| BWL-KCRP-600-60 | Fits 40" to 50" Diameter |