

SINCE 1925

KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS

Steam Component Bags with Drawstring

Description



- Designed for steam sterilization
- Convenient Drawstring Closure eliminates the need for autoclave tape or the validation and use of a heat sealer
- Fast and consistent for equipment preparation area
- Strong, sealed seams eliminates sewing threads and greatly reduce particulates
- Resists tearing and punctures
- Manufactured and packaged in a certified ISO 6 Cleanroom
- Made in the USA
- Certificates of Conformance sent with shipment

Tyvek® is a trademark of affiliates of DuPont de Nemours, Inc.

PROPERTY	UNITS	TEST METHOD	TYPICAL VALUES 1422A TYVEK®
Seal Strength	N/mm	F-88	>.78
Basis Weight	g/m ²	EN ISO 536	41
Total Thickness	μ	EN 20534	130
Air Permeability	s	ISO 5636-5	20
Water Vapor Permeability	G/M ² /24 HR	ASTM 398	1700
Hydrostatic Head	cm h ₂ o	ISO 811	120
Trapezoidal Tear MD/CD	N	ISO 9073-4	26.1/30.6
Linting Propensity	Rating	BS 6909 (Shirley Method 21)	Pass Rating: Excellent
Burst Resistance	kPa	ISO 2960 (50 cm ²)	108
Puncture Resistance	N	EN 863	10.8
Melting Point	c°		135

SINCE 1925



SAFETY ↔ CLEANROOM PRODUCTS

Size W x L	Keystone Part No.
1" x 1"	60101TWTD
1" x 4"	60105TWTD
1" x 5"	60106TWTD
1" x 6"	60107TWTD
1" x 7"	60108TWTD
1" x 8"	60109TWTD
1" x 11"	60112TWTD
1" x 17"	60117TWTD
1" x 20"	60120TWTD
1" x 24"	60124TWD
2" x 2"	60202TWTD
2" x 3"	60204TWTD
2" x 4"	60205TWTD
2" x 5"	60206TWTD
2" x 6"	60207TWTD
2" x 7"	60208TWTD
2" x 8"	60209TWTD
2" x 9"	60210TWTD

Size W x L	Keystone Part No.
2" x 10"	60211TWTD
2" x 11"	60212TWTD
2" x 12"	60212TWTD
2" x 13"	60213TWTD
2.5" x 3"	60406TWTD
3" x 5"	60306TWTD
3" x 6"	60307TWTD
3" x 8"	60709TWTD
3" x 9"	60310TWTD
3" x 10"	60311TWTD
3" x 12"	60313TWTD
3" x 14"	60314TWTD
4" x 3"	60404TWTD
4" x 4"	60405TWTD
4" x 7"	60407TWTD
4" x 9"	60409TWTD
4" x 10"	60411TWTD
4" x 12"	60412TWTD

SINCE 1925



SAFETY ↔ CLEANROOM PRODUCTS

Size W x L	Keystone Part No.
4" x 18"	60418TWTD
4" x 25"	60425TWTD
4" x 26"	60426TWTD
4" x 28"	60428TWTD
4" x 30"	60430TWTD
4" x 43"	6043TWTD
4" x 48"	60448TWTD
5" x 6"	60507TWTD
5" x 8"	60509TWTD
5" x 9"	60510TWTD
5" x 10"	60511TWTD
5" x 11"	60512TWTD
5" x 12"	60513TWTD
5" x 15"	60515TWTD
6" x 6"	60607TWTD
6" x 8"	60609TWTD
6" x 9"	60610TWTD
6" x 10"	60611TWTD

Size W x L	Keystone Part No.
6" x 12"	6012TWTD
6" x 15"	60615TWTD
6" x 16"	60616TWTD
6" x 17"	6017TWTD
6" x 23"	60623TWTD
6" x 33"	60633TWTD
6" x 36"	60636TWTD
6" x 38"	60638TWTD
6" x 43"	60643TWTD
6" x 62"	60662TWTD
7" x 7"	60708TWTD
7" x 9"	60710TWTD
7" x 10"	60711TWTD
7" x 13"	60713TWTD
7" x 10"	60710TWTD
7" x 17"	60717TWTD
7" x 27"	60727TWTD
8" x 8"	60809TWTD

SINCE 1925



SAFETY & CLEANROOM PRODUCTS

Size W x L	Keystone Part No.
8" x 11"	60812TWTD
8" x 12"	60813TWTD
8" x 15"	60815TWTD
8" x 17"	60817TWTD
8" x 22"	60822TWTD
8" x 24"	60824TWTD
8" x 31"	60831TWTD
8" x 43"	60843TWTD
8" x 56"	60856TWTD
8" x 65"	60865TWTD
8" x 72"	60872TWTD
9.5" x 12"	6000TWTD
10" x 15"	61015TWTD
10" x 20"	61020TWTD
10" x 24"	61024TWTD
10" x 27"	61027TWTD
10" x 30"	61030TWTD
10" x 37"	61037TWTD

Size W x L	Keystone Part No.
10" x 43"	61043TWTD
10" x 57"	61057TWTD
10" x 62"	61062TWTD
12" x 12"	61212TWTD
12" x 16"	61216TWTD
12" x 18"	61218TWTD
12" x 20"	61220TWTD
12" x 24"	61224TWTD
12" x 27"	61227TWTD
12" x 36"	61236TWTD
12" x 48"	61248TWTD
14" x 16"	61416TWTD
14" x 18"	61418TWTD
14" x 20"	61420TWTD
14" x 22"	61422TWTD
14" x 31"	601422TWTD
15" x 19"	61519TWTD
15.5" x 19.5"	6155195TWTD

SINCE 1925



SAFETY ↔ CLEANROOM PRODUCTS

Size W x L	Keystone Part No.
15" x 26"	61526TWTD
15" x 30"	61530TWTD
15" x 44"	61544TWTD
16" x 16"	61616TWTD
16" x 18"	61618TWTD
16" x 24"	6002TWTD
16" x 27"	61627TWTD
16" x 30"	61630TWTD
16" x 32"	61632TWTD
16" x 38"	61638TWTD
16" x 42"	61642TWTD
18" x 19"	61819TWTD
18" x 24"	61824TWTD
18" x 30"	6003TWTD
18" x 36"	61836TWTD
18" x 48"	61848TWTD
18" x 54"	61854TWTD
18" x 55"	61855TWTD

SINCE 1925



SAFETY ↔ CLEANROOM PRODUCTS

Size W x L	Keystone Part No.
18" x 70"	61870TWTD
20" x 20"	62020TWTD
20" x 22"	62022TWTD
20" x 24"	62024TWTD
20" x 26"	62026TWTD
20" x 30"	62030TWTD
20" x 33"	62033TWTD
20" x 37"	62037TWTD
20" x 40"	6004TWTD
22" x 18"	62218TWTD
22" x 22"	62222TWTD
22" x 25"	62225TWTD
22" x 30"	62230TWTD
23" x 25"	62325TWTD
24" x 24"	6008TWTD
24" x 30"	62430TWTD
24" x 36"	6005TWTD
24" x 40"	62440TWTD

SINCE 1925



SAFETY ⇄ CLEANROOM PRODUCTS

Size W x L	Keystone Part No.
24" x 42"	62442TWTD
24" x 48"	62448TWTD
24" x 55"	62455TWTD
24" x 62"	62462TWTD
26" x 32"	62632TWTD
26" x 34"	62634TWTD
26" x 42"	62642TWTD
28" x 18"	62818TWTD
28" x 28"	62828TWTD
28" x 34"	62834TWTD
28" x 38"	62838TWTD
28" x 44"	62844TWTD
30" x 30"	63030TWTD
30" x 36"	63036TWTD
30" x 40"	63040TWTD
30" x 45"	63045TWTD
30" x 53"	63053TWTD
30" x 60"	63060TWTD

Size W x L	Keystone Part No.
32" x 46"	63246TWTD
32" x 53"	63253TWTD
32" x 58"	63258TWTD
32" x 63"	63263TWTD
34" x 34"	63434TWTD
34" x 37"	63437TWTD
34" x 41"	63441TWTD
34" x 48"	63448TWTD
36" x 36"	63636TWTD
36" x 38"	63638TWTD
36" x 43"	63643TWTD
36" x 46"	63646TWTD
36" x 48"	63648TWTD
36" x 64"	63664TWTD
40" x 40"	604040TWTD
40" x 47"	604047TWTD
40" x 49"	64049TWTD
40" x 52"	64052TWTD

SINCE 1925

KEYSTONE

SAFETY ⇄ CLEANROOM PRODUCTS

Size W x L	Keystone Part No.
40" x 60"	64060TWTD
45" x 45"	64545TWTD
48" x 48"	64848TWTD
48" x 60"	64860TWTD
48" x 64"	64864TWTD

SINCE 1925

KEYSTONE

SAFETY ⇄ CLEANROOM PRODUCTS