

SINCE 1925

KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS

Steam Sterilization Wrapper Made with 1422A Tyvek®



Description

- Designed for steam sterilization
- Recommended for use in conjunction with a secondary sterilization pouch made with Tyvek®
- Hydrophobic
- Excellent microbial barrier
- Resists tearing and punctures
- Lint-Free
- Manufactured and packaged in ISO certified cleanrooms
- Certificates of Conformance sent with all shipments
- Made in USA

Tyvek® is a trademark of affiliates of DuPont de Nemours, Inc.

PROPERTY	UNITS	TEST METHOD	TYPICAL VALUES 1422A TYVEK®
Basis Weight	g/m ²	EN ISO 536	41
Total Thickness		EN 20534	130
Air Permeability	s	ISO 5636-5	20
Water Vapor Permeability	G/M ² /24 HR	ASTM 398	1700
Hydrostatic Head	cm h ₂ O	ISO 811	120
Trapezoidal Tear MD/CD	N	ISO 9073-4	26.1/30.6
Linting Propensity	Rating	BS 6909 (Shirley Method 21)	Pass Rating: Excellent
Burst Resistance	kPa	ISO 2960 (50 cm ²)	108
Puncture Resistance	N	En 863	10.8
Melting Point	C		135

SINCE 1925

KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS

Dimensions	Keystone Part Number
3 x 3 in.	T2ASHEET0303
4 x 4 in.	T2ASHEET0404
6 x 6 in.	T2ASHEET0606
7.5 x 7.5 in.	T2ASHEET07-507-5
9 x 9 in.	T2ASHEET0909
9 x 11 in.	T2ASHEET0911
12 x 12 in.	T2ASHEET1212
15 x 15 in.	T2ASHEET1515
18 x 18 in.	T2ASHEET1818
20 x 20 in.	T2ASHEET2020
24 x 24 in.	T2ASHEET2424
30 x 30 in.	T2ASHEET3030
40 x 40 in.	T2ASHEET4040
40 x 80 in.	T2ASHEET4080
60 x 60 in.	T2ASHEET6060

SINCE 1925



SAFETY ↔ CLEANROOM PRODUCTS

Dimensions	Keystone Part Number
60 x 140 in.	T2ASHEET60140
72 x 72 in.	T2ASHEET7272
84 x 94 in.	T2ASHEET8494

SINCE 1925



SAFETY ↔ CLEANROOM PRODUCTS