


SINCE 1925

KEYSTONE

SAFETY ⇄ CLEANROOM PRODUCTS

Stopper Track Cover made with 1073B Tyvek®



Description

- Tyvek® construction for excellent strength and exceptional microbial barrier
- Designed for steam sterilization
- Bendable, flexible
- Recommended for sterilization of stopper track
- Unlike elastic versions, designed to remain in place during stopper bowl and track set up
- Manufactured and packaged in ISO certified cleanroom
- Improves the efficiency of the parts preparation department
- Certificates of Conformance sent with each shipment
- Made in USA

Tyvek® is a trademark of affiliates of DuPont de Nemours, Inc.

SINCE 1925

KEYSTONE

SAFETY ✦ CLEANROOM PRODUCTS

PROPERTY	UNITS	TEST METHOD	TYPICAL VALUES 1073B TYVEK®	TYPICAL VALUES PET/POLYOLEFIN FILM
Seal Strength	Lbs/inchwidth	F-88	>1.5	>1.5
Basis Weight	oz/yd2	ASTM D3776 EN ISO 536	2.2	N/A
Total Thickness	Mils	D 1777-64	7	2.6
Microbial Barrier	Log Reduction Value (LRV)	ASTM F 1608	5.2	N/A
Bendston Air Permeability	mL/min	ISO 5636-3	572	N/A
Moisture Vapor Transmission Rate	g/m2/24 hr	TAPPI T523	1615	N/A
Hydrostatic Head	Inches H2O	AATCC TM 127 EN20811	58	N/A
Tensile Strength MD	lb/in.	ASTM D5035 EN ISO 1924-2	44	N/A
Tensile Strength CD	lb/in.	ASTM D5035 EN ISO 1924-2	45	N/A
Elongation MD	%	ASTM D5035 EN ISO 1924-2	20	N/A
Elongation CD	%	ASTM D5035 EN ISO 1924-2	24	N/A
Elmendorf Tear, MD	lbs	ASTM 1424 EN 21974	0.7	N/A
Elmendorf Tear, CD	lbs	ASTM 1424 EN 21974	0.8	N/A
Mullen Burst	psi	ASTM D774 ISO 2758	176	N/A
Spencer Puncture	in.-lb f/in.2	ASTM 3420	50	N/A
Opacity	%	TAPPI T425 ISO 2471	91	6

Size W x L	Keystone Part No.
16" x 8"	MGTVKSHEET1608WT4
22" x 16"	MGTVKSHEET221610WT4-2PLY
28" x 10"	MGTVKSHEET2810WT4-2PLY
34" x 14"	MGTVKSHEET3414WT8
36" x 17"	MGTVKSHEET3617WT8
42" x 17"	MGTVKSHEET4217WT8

SINCE 1925

KEYSTONE

SAFETY ✦ CLEANROOM PRODUCTS