

SINCE 1925

KEYSTONE

SAFETY ⇄ CLEANROOM PRODUCTS

Stopper Track Cover made with 1421B Tyvek®



Description

- Uncoated Tyvek® for maximum breathability during sterilization without the chance of residue deposits on equipment, post sterilization
- Designed for steam sterilization
- Bendable, flexible
- Recommended for sterilization of stopper track
- Unlike elastic versions, designed to remain in place during stopper bowl and track set up
- Manufactured and packaged in ISO certified cleanroom
- Improves the efficiency of the parts preparation department
- Certificates of Conformance sent with each shipment
- Made in USA

Tyvek® is a trademark of affiliates of DuPont de Nemours, Inc.

Property	Comparable standard test method	Units	Nominal value (range) Tyvek® 1421B
Preliminary product specifications (English units)			
Basis weight	ASTM D3776 ¹	oz/yd ²	1.22 [1.16–1.26]
Delamination	ASTMD2724 ²	lb _f	0.08 [>0.03]
Gurley Hill porosity	TAPPI T460 ³	sec/100 cc	20 [8–36]
Preliminary product specifications (metric units)			
Basis weight	ASTM D3776 ¹	g/m ²	41 [39–43]
Delamination	ASTMD2724 ²	N	0.36 [>0.13]
Gurley Hill porosity	TAPPI T460 ³	sec/100 cc	20 [8–36]
Preliminary miscellaneous properties (English units)			
Thickness	EN ISO 534 ¹	mils	6 (0.15mm)
Tensile strength, MD	EN ISO 1924-2 ²	lb _f /in	7.6
Tensile strength, CD	EN ISO 1924-2 ²	lb _f /in	6.4
Elmendorf tear, MD	ASTM D1424	lb _f	2.4
Elmendorf tear, CD	ASTM D1424	lb _f	1.9
Hydrostatic head	AATCC TM 127 ³	in H ₂ O	44
Elongation, MD	EN ISO 1924-2 ²	%	10
Elongation, CD	EN ISO 1924-2 ²	%	15
Puncture resistance	EN 863 ⁴	lb _f	2.7
Mullen burst	ISO 2758	psi	53
Bendsten air permeability	ISO 5636-3 ⁵	mL/min	741
Microbial barrier	ASTM F1608	LRV	~4

MD = machine direction; CD = cross direction

1. Surface 2 cm², pressure 50kPa.

2. Modified for speed, sample width 2.54 cm (1 inch).

3. Rate of use: 60 cm H₂O/min.

4. Modified for conditioning.

5. Electronic device.

SINCE 1925

KEYSTONE

SAFETY ✦ CLEANROOM PRODUCTS

Size W x L	Keystone Part No.
4" x 10"	70TK0410WT2-2PLY
35" x 17"	70TK3517WT8
36" x 17"	70TK3617WT8

SINCE 1925

KEYSTONE

SAFETY ✦ CLEANROOM PRODUCTS