

SINCE 1925

KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS

Stopper Bowl Cover, Made with Tyvek®

Save time by using the stopper bowl cover instead of wrapping with sheets of material and autoclave tape. When used in conjunction with a secondary Tyvek® Bag, it creates an effective wrapping system. After the outer bag is removed, this stopper bowl cover should remain in place, providing protection from microbial and particulate contamination while the remaining set-up activities are performed. Tyvek® is a trademark of affiliates of DuPont de Nemours, Inc.



Typical Application & Characteristics:

- Steam sterilization indicator
- Autoclavable at 121°C-125°C
- Hydrophobic
- Resists tearing and punctures
- Lint-free
- Manufactured and packaged in ISO certified cleanrooms
- Conveniently packaged: 1 per heat-sealed and labeled poly bag
- Certificates of Conformance sent with all shipments
- Made in USA

PROPERTY	UNITS	TEST METHOD	TYPICAL VALUES 1422A TYVEK®
Linting Propensity	Rating	BS 6909 (Shirley Method 21)	Pass Rating: Excellent
Basis Weight	g/m ²	EN ISO 536	41
Total Thickness	μ	EN 20534	130
Air Permeability	s	ISO 5636-5	20
Water Vapor Permeability	G/M/24 HR	ASTM 398	1700
Hydrostatic Head	cm H ₂ O	ISO 811	120
Trapezoidal Tear MD	N	ISO 9073-4	26.1
Trapezoidal Tear CD	N	ISO 9073-4	30.6
Burst Resistance	kPa	ISO 2960 (50 cm ²)	108
Puncture Resistance	N	EN 863	10.8
Melting Point	C [°]		135

SINCE 1925

KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS



DIMENSIONS	KEYSTONE PART NO.
Fits 11" to 16" Diameter	BWLT2A20
Fits 12" to 20" Diameter	BWLT2A24
Fits 14" to 22" Diameter	BWLT2A26
Fits 20" to 26" Diameter	BWLT2A30
Fits 26" to 30" Diameter	BWLT2A36
Fits 28" to 35" Diameter	BWLT2A42
Fits 30" to 38" Diameter	BWLT2A46
Fits 30" to 40" Diameter	BWLT2A54
Fits 35" to 45" Diameter	BWLT2A54
Fits 40" to 50" Diameter	BWLT2A60
Fits 50" to 60" Diameter	BWLT2A72