

SINCE 1925

KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS

UV LIGHT BLOCKING TYVEK® FC COVERS



UV Light Blocking Tyvek® FC Equipment Covers are designed to keep cleaned pharmaceutical equipment and supplies clean and protected from external sources of contamination while also providing in-process protection to light sensitive biotech products.

TYPICAL APPLICATION & CHARACTERISTICS:

- Each cover is individually poly bagged and labeled for easy identification
- Available with, or without elastic with, or without zippers, or snaps
- UV Light Blocking Tyvek® FC fabric protects light sensitive Biotech products
- Hydrophobic
- Resists tearing and punctures
- Anti-static
- Manufactured and packaged in ISO certified Cleanroom
- Certificates of Conformance sent with each shipment
- Made in USA

MATERIAL SPECIFICATIONS:

PROPERTY	UNITS	TEST METHOD	TYPICAL VALUES TYVEK® 400 FC
Basis Weight	oz/yd ²	ASTM D3776	1.6 oz/yd ²
Thickness	mils	ASTM D1117	7.4 mils
Breaking Strength-Grab (MD)	f	ASTM D5034	20 lbf
Breaking Strength-Grab (CD)	f	ASTM D5034	26 lbf
Hydrostatic Head	in H ₂ O	AATCC 127	4.2 in H ₂ O
Tear Resistance MD	f	ASTM D5733	5lbf
Tear Resistance CD	f	ASTM D5733	7lbf
Flammability		16 CFR 1610	Class 1

SINCE 1925

KEYSTONE

SAFETY ↔ CLEANROOM PRODUCTS

UV LIGHT BLOCKING TYVEK® FC COVERS



KEYSTONE PART NUMBER	DIMENSIONS
COV-GTVK-141008	14" x 10" x 8"
COV-GTVK-262010	26" x 20" x 10"
COV-GTVK-301616	30" x 16" x 16"
COV-GTVK-322212	30" x 22" x 12"
COV-GTVK-302456	30" x 24" x 56"
COV-GTVK-322824	32" x 28" x 24"
COV-GTVK-362660	36" x 26" x 60"
COV-GTVK-423460	42" x 34" x 60"
COV-GTVK-442430	44" x 24" x 30"
COV-GTVK-443632	44" x 36" x 32"
COV-GTVK-444460	44" x 44" x 60"
COV-GTVK-482660	48" x 26" x 60"
COV-GTVK-482872	48" x 28" x 72"
COV-GTVK-484860	48" x 48" x 60"
COV-GTVK-502474	50" x 24" x 74"
COV-GTVK-502872	50" x 28" x 72"
COV-GTVK-505050	50" x 50" x 50"
COV-GTVK-545486	54" x 54" x 86"
COV-GTVK-563660	56" x 36" x 60"
COV-GTVK-602460	60" x 24" x 60"
COV-GTVK-602854	60" x 28" x 54"
COV-GTVK-603072	60" x 30" x 72"
COV-GTVK-604444	60" x 44" x 44"
COV-GTVK-605060	60" x 50" x 60"
COV-GTVK-642460	64" x 24" x 60"
COV-GTVK-642478	64" x 24" x 78"
COV-GTVK-643068	64" x 30" x 68"
COV-GTVK-654058	65" x 40" x 58"
COV-GTVK-682238	68" x 22" x 38"
COV-GTVK-742686	74" x 26" x 86"
COV-GTVK-762672	76" x 26" x 72"

